


# 96-Well LOFT Plate<sup>®</sup>

for kinetic endotoxin testing



Save 85%  
reagents



Prevent  
contamination



Eliminate  
evaporation



Improve  
accuracy



Meet USP, EP, and JP  
standards


# Easy Operation Steps



Add sample




Centrifuge



Readout

## Innovative LOFT Plate®

It is an endotoxin plate at another level.



1. Save 85% LAL and endotoxin reagent consumption
2. Skim air bubbles
3. Improve readout accuracy by a fixed optical path
4. Assist mixing by centrifugation
5. Eliminate evaporation and environmental contamination by dual optical windows
6. Compatible with existing endotoxin test system